



<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Date</b>
ScImage, Inc.	PicomEnterprise	2.0	January 20, 2006
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>		<b>Options Implemented</b>
Scheduled Workflow	Image Manager/Image Archive		
	Performed Procedure Step Manager		
Patient Information Reconciliation	Image Manager/Image Archive		
Cardiac Catheterization Workflow	Image Manager/Image Archive		Cardiac CATH Intermittently Connected Modality
	Performed Procedure Step Manager		
ECHO Workflow	Image Manager/Image Archive		Echocardiography Intermittently Connected Modality
	Performed Procedure Step Manager		
Consistent Time	Time Client		SNTP
<b>Internet address for vendor's IHE information:</b> <a href="http://www.scimage.com/conformance.asp">www.scimage.com/conformance.asp</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	<a href="http://www.scimage.com/hl7.asp">www.scimage.com/hl7.asp</a>		
<b>DICOM</b>	<a href="http://www.scimage.com/scimagedcs.asp">www.scimage.com/scimagedcs.asp</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.rsna.org/IHE">www.rsna.org/IHE</a>	In Europe: <a href="http://www.ihe-europe.org">www.ihe-europe.org</a>	In Japan: <a href="http://www.jira-net.or.jp/ihe-j">www.jira-net.or.jp/ihe-j</a>	